

# Impact of Panel Size on Molecular Residual Disease (MRD) Assay Performance

C.J. Battey, PhD\*; Ashley Acevedo, PhD\*; Matt LaBella, PhD\*; Sangita Ganesh, PhD; Elise Buser; Nafei Xu; Ravi Patel, PhD; Kiefer Haug; Kyle Trettin, PhD; Jeff Jasper, PhD; Jan Risinger; Genevieve Gould, PhD; Dale Muzzey, PhD

\*Authors contributed equally

Myriad Genetics, Inc. Salt Lake City, UT, USA



## Background

- Sequencing cell-free DNA (cfDNA) is a promising method for monitoring cancer treatment response and detecting recurrence. However, sensitivity at the low tumor fractions typical of early-stage, post-treatment, and early-recurrent tumors is limited by the small number of variants targeted in commercially available assays.
- We developed a tumor-informed MRD assay using tumor-normal whole-genome sequencing (WGS) followed by interrogation of select somatic variants in cfDNA. The method was tested on a cohort of 30 patient samples, targeting 1,000 somatic single-nucleotide variants (SNV) per sample.

## WGS Panel Design

- 20-50x tumor-normal WGS data was processed with a bioinformatic pipeline including somatic SNV calling, tumor purity estimation, and copy number anomaly calling.
- 1,000 target SNVs from each sample were selected to create a sample-specific panel for hybridization capture and deep sequencing.
- When fewer than 1,000 high-confidence targets were identified we backfilled with lower-confidence variants.
- We evaluated somatic calling performance by deep targeted sequencing of formalin-fixed paraffin-embedded (FFPE) tumor tissue (Figure 1).

## Circulating Tumor DNA Detection & Quantification

- Somatic targets were enriched and sequenced to high depth (>500x after UMI deduplication) in FFPE tumor tissue, patient normal DNA, and pre-treatment patient cfDNA.
- We developed a maximum likelihood statistical model to quantify tumor fraction (the genome-equivalents proportion of cfDNA derived from tumor tissue) and call MRD positive or negative status in patient cfDNA.
- We then created serial dilutions of patient cfDNA with non-patient cfDNA at a set of predetermined tumor fractions for 17 samples with sufficient cfDNA mass and tumor fraction.
- Bootstrap resampling of target sites was used to estimate sensitivity across the tumor fraction spectrum given varying target site counts (Figure 2).
- Tumor-fraction quantification accuracy was estimated by comparing model estimates to the expected dilution tumor fraction (Figure 3).

## Conclusions

- WGS-driven panel design allows targeting up to 1,000 high-confidence somatic variants in an MRD assay across diverse cancer types and stages.
- Increasing the MRD panel size in turn leads to higher sensitivity at low tumor fractions and robust tumor-fraction quantification in patient cfDNA.
- High-sensitivity MRD has the potential to enable earlier recurrence detection and gives researchers and clinicians new tools for monitoring patient treatment responses.

## Results

Figure 1. WGS Panel Design Variant Counts and Target Confirmation Rate




Figure 2. Assay Sensitivity as a Function of Tumor Fraction and Target Variant Count in 17 cfDNA Dilution Series

Models are tuned for 100% specificity across 136 negative controls.




Figure 3. Tumor-Fraction Quantification in 17 cfDNA Dilution Series

Tumor-fraction quantification was linear down to 0.001%, with a mean absolute error of 0.0016%.

